




1-1-1

r AyII Migraine






rAyII

Migraine (Spring)



<i>rAyI/0:6</i>	<i>rTyII/0:c</i>
	

rAyII Migraine (Summer)

<i>rAyI/0:6</i>	<i>rAyIII/0:6</i>	<i>rTxIII/1:c</i>
		

r AyII



Migraine (Fall)

<i>lAxI/1:(6)</i>	<i>lAxII/2:c</i>
	

r AyII Migraine (Winter)




<i>lAxII/2:c</i>	<i>rAyII/0:c</i>
	

rAyII Migraine (Dull Pain)

<i>rTxIII/1:c</i>	<i>ЛxI/1:(c)</i>
	

rAyII

Migraine (Muscle Pain)

<i>rAyI/0:6+c</i>	<i>rAyIII/0:6+c</i>	<i>rAyII/0:3+c</i>
		

r AyII

Migraine

Spring

rAyI/0:6

C point not used, F point "6" using No. 3 head, stimulate proximally.

rTyII/0:c

C point not used, F point "c" using No. 2 head, stimulate distally.

Summer

rAyI/0:6

C point not used, F point "6" using No. 3 head, stimulate proximally.

rAyIII/0:6

C point not used, F point "6" using No. 3 head, stimulate proximally.

rTxIII/1:c

C point "1" using No. 1 head, F point "c" using No. 2 head, stimulate proximally.

Fall

ITxI/1:(6)

C point "1" using No. 3 head, F point "6" using No. 3 head, stimulate distally.

lAxII/2:c

C point "2" using No. 3 head, F point "c" using No. 2 head, stimulate distally.

Winter

lAxII/2:c

C point "2" using No. 3 head, F point "c" using No. 2 head, stimulate distally.

rAyII/0:c

C point not used, F point "c" using No. 2 head, stimulate proximally.

Dull Pain

rTxIII/1:c

C point "1" using No. 1 head, F point "c" using No. 2 head, stimulate proximally.

ITxI/1:(c)

C point "1" using No. 3 head, F point "c" using No. 2 head, stimulate distally.

Muscle Pain

rAyI/0:6+c

C point not used, F point "6" using No. 3 head, "c" using No. 2 head, stimulate proximally.

rAyIII/0:6+c

C point not used, F point "6" using No. 3 head, "c" using No. 2 head, stimulate proximally.

rAyII/0:3+c

C point not used, F point "3" using No. 3 head, "c" using No. 2 head, stimulate proximally.

rAyII

偏頭痛治療方法

季節	痠痛、沈重感、肌肉痛
<p>春：</p> <ol style="list-style-type: none"> rAyI/0:6 不使用C按點，F按點「6」使用3號壓棒、向近端方向刺激。 rTyII/0:c 不使用C按點，F按點「c」使用2號壓棒、向遠端方向刺激。 <p>夏：</p> <ol style="list-style-type: none"> rAyI/0:6 不使用C按點，F按點「6」使用3號壓棒、向近端方向刺激。 rAyIII/0:6 不使用C按點，F按點「6」使用3號壓棒、向近端方向刺激。 rTxIII/1:c C按點「1」使用1號壓棒，F按點「c」使用2號壓棒、向近端方向刺激。 <p>秋：</p> <ol style="list-style-type: none"> ITxI/1:(6) C按點「1」使用3號壓棒，F按點「6」使用3號壓棒、向遠端方向刺激。 IAxII/2:c C按點「2」使用3號壓棒，F按點「c」使用2號壓棒、向遠端方向刺激。 <p>冬：</p> <ol style="list-style-type: none"> IAxII/2:c C按點「2」使用3號壓棒，F按點「c」使用2號壓棒、向遠端方向刺激。 rAyII/0:c 不使用C按點，F按點「c」使用2號壓棒、向近端方向刺激。 	<p>痠痛、沉重感：</p> <ol style="list-style-type: none"> rTxIII/1:c C按點「1」使用1號壓棒，F按點「c」使用2號壓棒、向近端方向刺激。 ITxI/1:(c) C按點「1」使用3號壓棒，F按點「c」使用2號壓棒、向遠端方向刺激。 <p>肌肉痛：</p> <ol style="list-style-type: none"> rAyI/0:6+c 不使用C按點，F按點「6」使用3號壓棒、「c」使用2號壓棒、向近端方向刺激。 rAyIII/0:6+c 不使用C按點，F按點「6」使用3號壓棒、「c」使用2號壓棒、向近端方向刺激。 rAyII/0:3+ c 不使用C按點，F按點「3」使用3號壓棒、「c」使用2號壓棒、向近端方向刺激。